


Low and medium energy Shielded hot cell

Technetium shielded isolator

- 1 shielded compartment for 2 generators with independent lifts
- 1 shielded compartment for dose calibrator
- 1 shielded compartment for waste container
- 2 glove ports
- 1 transfer chamber



Available with preparation robot

- Doubling operator's preparation capacity
- Pharmaceutical quality and traceability
- Reduced radiation exposure

For further information:
sales@radioprotech.com - +33 770 199 101
More products: www.radioprotech.com



High energy Shielded hot cell



Gallium shielded isolator

- 1 shielded compartment for 2 Ge / Ga generators
- 1 shielded compartment for dose calibrator
- 1 shielded compartment for waste container
- 1 transfer chamber
- 2 glove ports

Research shielded isolator

- 1 shielded compartment for FDG container
- 1 shielded compartment for 2 Ge / Ga compartment
- 1 shielded compartment for waste container
- 1 transfer chamber
- 4 glove ports

