

Shielded trolley

Iodine-131 capsules transport trolley

To collect, store and transport radioactive waste

- Shielding: 15 to 40 mm lead




Tc99 generator transport trolley

To transport Tc99 generators including uncontrolled areas

- Shielding: 15 to 40 mm lead

All our trolleys can be motorized




For further information:
sales@radioprotech.com - +33 770 199 101
More products: www.radioprotech.com

Waste management

Shielded waste container

To collect and store solid radioactive waste

- Shielding: 5, 10 or 20 mm lead
- Capacity: 30, 60 and 80 liters
- For low, medium and high energy



Shielded sharp container

For main clinical waste containers

- Size standard or XL from 1,5 to 3,5 liters
- Shielding: 3, 5, 10 or 15 mm lead

For further information:

sales@radioprotech.com - +33 770 199 101

More products: www.radioprotech.com