

Injection shield


For high energy

- Shielded by 30 mm lead
- Large window
- Height adjustable shelf
- PET / MR compatible




Options:

Infusion bag holder



Tele-injector



For further information:
sales@radioprotech.com - +33 770 199 101
More products: www.radioprotech.com






Tabletop shield

Gamma emitters

- From 5 to 30 mm lead
- From 20 to 60 mm lead glass



Mobile tabletop



Gamma and beta emitters

- 10 mm PMMA
- 5 mm lead
- 20 mm lead glass



Beta emitters

- 10 mm PMMA

For further information:
sales@radioprotech.com - +33 770 199 101
More products: www.radioprotech.com

