


Shielded biosafety cabinet

For radiopharmaceuticals preparation or cells labeling




- Class II biosafety cabinet
- Sliding lead glass window
- Several levels of equipment
- Shielding: 5 or 10 mm lead

For further information:
sales@radioprotech.com - +33 770 199 101
More products: www.radioprotech.com



Shielded biosafety cabinet





Generators compartment

- Shielded with 50 mm lead
- Compatible with all generators and brands
- 2 independents lifts

Dose calibrator compartment

- Shielded with 10 mm lead
- Compatible with all dose calibrators and brands foot-operated elevator



Waste compartment

- Shielded with 10 mm lead
- 2 bins containers
- 2 filling holes

For further information:

sales@radioprotech.com - +33 770 199 101

More products: www.radioprotech.com

